Little Doctor

LD-60, LD-61, LD-70, LD-70NR, LD-71, LD-71A, LD-80, LD-100

ENG BLOOD PRESSURE MONITOR. Instruction Manual.
POL CIŚNIENIOMIERNYK MECHANICZNY LD DO POMIARU CIŚNIENIA TĘTNICZEGO KRWI. Instrukcja obsługi.
HUN LD MŰSZER AZ ARTERIÁS VÉRNYOMÁS MÉRÉSÉRE. Használati útmutató.
ROU MONITOR DE TENSIUNE ARTERIALĂ LD. Manual de instrucțiuni.
BGR УСТРОЙСТВО ЗА ИЗМЕРВАНЕ НА КРЪВНОТО НАЛЯГАНЕ LD. Ръководство за експлоатация.

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ENG

PARTS AND COMPONENTS (fig.1)

№	Name of part	Model	Quantity							
			LD-60	LD-61	LD-70	LD-70NR	LD-71	LD-71A	LD-80	LD-100
1	Manometer	Manometer LD-S013	1	1	1	1	1	1	1	–
2	Manometer	Manometer LD-S048	–	–	–	–	–	–	–	1
3	Adult's cuff with ring	Cuff NZAR Pneumatic chamber LD-S02A	–	–	1	–	1	1	–	1
	Adult's cuff	Cuff N2A Pneumatic chamber LD-S02A	–	1	–	1	–	–	–	–
4	Big adult's cuff with ring	Cuff NZLR Pneumatic chamber LD-S02L	1	–	–	–	–	–	–	–
	Child's cuff	Cuff CZC Pneumatic chamber LD-S02C	–	–	–	–	–	–	1	–
5	Cuff for newborn infants	Cuff CZM Pneumatic chamber LD-S02N	–	–	–	–	–	–	1	–
	Cuff for infants	Cuff CZI Pneumatic chamber LD-S02I	–	–	–	–	–	–	1	–
6	Stethoscope with built-in head into cuff	LD Prof-Plus (with LD-S025a head)	1	–	–	–	–	1	–	–
	Stethoscope	LD Prof-Plus (with LD-S025 head)	–	1	–	–	–	1	–	–
7	Air pump assembly	Bulb LD-S014 Air valve LD-S015 Non-return valve LD-S016	1	1	1	1	1	1	1	1
	Storage Case	Storage Case LD-S059	1	1	1	1	1	1	1	–
8	Instruction Manual + Warranty Card		–	1	1	1	1	1	1	1

GENERAL INFORMATION

This Instruction Manual is designed to assist the user with safe and effective operation of the Device for measurement of blood pressure (hereinafter – the "Device"). Use this Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual, in particular «RECOMMENDATIONS ON CORRECT MEASUREMENT».

INDICATIONS FOR USE

The device is designed to measure arterial pressure of human according to Korotkov method. The device is recommended for use under conditions of clinics and hospitals, as well as in household conditions as supplement to medical observation. Pressure is measured by means of auscultation of Korotkov's tones by stethoscope and by taking readings on manometer.

RECOMMENDATIONS ON CORRECT MEASUREMENTS

- Do not use the device without preliminary consultation with Your doctor, if You undergo treatment by hemodialysis or by anticoagulants, antithrombocytes or steroids. Use of device in these cases may cause internal hemorrhage.
- For correct measurement you should know that THE BLOOD PRESSURE IS SUBJECT TO SHARP VARIATIONS EVEN WITHIN THE SHORT TIME INTERVALS. The blood pressure depends on many factors. It is usually lower in summer and higher in winter. The blood pressure varies together with the atmospheric pressure, depends on physical loads, emotional excitement, stresses and dietary regime. Drugs, drinking alcohol and smoking produce significant effect. Even the very procedure of blood pressure measurement in a polyclinic sends the blood pressure high in many people, thus, the blood pressure measured at home often differs from the readings received in a polyclinic. As the blood pressure tends to rise at low temperatures, make measurements at an indoor temperature (approximately 20°C). If this Device stayed under a low temperature, keep it for at least 1 hour at an indoor temperature before use, otherwise the measurement result may be incorrect. During a day the difference in readings for healthy people may be 30-50 mmHg of systolic pressure and to 10 mmHg of diastolic pressure. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with blood pressure records. ONLY A CERTIFIED DOCTOR USING YOUR RECORDS IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VARIATIONS.
- At cardiovascular and some other diseases requiring blood pressure monitoring make measurements in the hours fixed by your attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. TAKING OF DRUGS AND THEIR DOSES SHOULD BE PRESCRIBED ONLY BY YOUR ATTENDING DOCTOR.
- KEEP QUIET DURING A MEASUREMENT TO OBTAIN THE ACCURATE VALUES OF YOUR BLOOD PRESSURE WITH THE ELECTRONIC DEVICE. Measure your blood pressure in the calm and comfortable conditions at the indoor temperature. No eating an hour before measurement; no smoking, taking tonic agents, alcohol 1-2 hours before measurement.
- The accuracy of blood pressure measurement depends on whether the cuff matches the size of your arm. THE CUFF SHOULD NOT BE TOO SMALL OR TOO LARGE.
- Wait 3 minutes between measurements for the blood to restore its circulation. However, the persons with prominent atherosclerosis due to considerable loss of vascular elasticity may need to increase the wait time between measurements (10-15 minutes). This also refers to the patients suffering for long from diabetes. For more accurate determination of blood pressure it is recommended to make a series of 3 consecutive measurements and to use the average value.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CUFF PREPARATION

- Insert the cuff end for about 5 cm into a metal ring as shown in the fig. 4.1 (if cuff is with ring).
- Apply the cuff to your left upper arm so that the air tube is directed to your palm (fig. 4.2). If the measurement is made on the left arm it is difficult, you may use your right arm. In this case remember that the readings may differ by 5-10 mmHg and even more.
- Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign "ARTERY" should be over the arm artery as shown in the fig. 4.3.
- Fix the cuff so that it fits tightly to the arm, but see that it is not overtight (fig. 4.4). Too tight or too free placement of the cuff may give inaccurate readings.
- On the fixed cuff the sign «INDEX» should point to the area «NORMAL» (fig. 4.5). It means that the cuff is chosen correctly and fits the size of your upper arm. If the sign points to the area marked «H» or to the left, the cuff is too small and the readings will be higher. If the sign points to the area marked «L» or to the right, the cuff is too large and the readings will be lower.
- If the arm has a conic form, the cuff should be put on with a spiral movement as shown in the fig. 4.6.
- If the rolled-up sleeve squeezes the arm interfering with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure (fig. 4.7).
- Place stethoscope head so that it was located in depression, somewhat higher than elbow bend (fig. 4.8). If you use the device with built-in stethoscope (for models LD-60 and LD-71A), to reach better auscultation of pulse on full hand, the cuff may be turned round the hand approximately by quarter of circle (60-90°) in such way that stethoscope head would be positioned on inner side of hand (more close to the body).

Pressure shall be measured in sitting or lying position of human.
IN SITTING POSITION WATCH THAT PART OF HAND WITH CUFF WOULD BE LOCATED AT THE LEVEL OF HEART, AND HAND WOULD BE FREELY LOCATED ON THE TABLE AND WOULD NOT MOVE.

CORRECT POSITION DURING MEASUREMENT (fig.5)

- Sit at a table so that during blood pressure measurement your hand rests on its surface. Be sure that the cuff is placed approximately at the level of your heart and that your arm lies freely on the table and does not move.
- You can measure the blood pressure lying on the back. Look at the ceiling, keep quiet and do not move during measurement. Be sure that the cuff is placed approximately at the level of your heart.

MEASUREMENT PROCEDURE

- Insert metal tubes of stethoscope into ears. Close air valve on bulb, having turned it clockwise. Squeezing bulb, pump up the cuff, listening to pulse by stethoscope. After You stopped to listen to pulse, pump up the cuff by 30 mm Hg more.
- Slowly opening air valve, turning it counterclockwise, bleed down pressure in the cuff. Watch that pressure in the cuff would fall down at speed 2-4 mm Hg per second. This is necessary to get accurate result.
- As soon as You hear weak pulse beats, memorize manometer reading. This is Your systolic (upper) arterial pressure.

- Pressure in the cuff is continuing falling down at the same speed (2-4 mm Hg per second). You are continuing to listen to pulse. Sounds which You hear will change. Unlike the first beats, they will become softer, resembling rustling. At the moment when You practically stop seize pulse, memorize manometer reading. This is Your diastolic (lower) arterial pressure.

CARE, STORAGE, REPAIR AND DISPOSAL

- Keep this Device from exposure to higher humidity, direct sunlight, shocks.
- Do not keep and use this Device near heating installations and open fire.
- Keep the Device clean and protect it from dust.
- Contact of device with aggressive solutions is not allowed.
- Protect the arm cuff and air tubes from contacting on sharp things.
- This Device does not contain special controls to adjust the measurement accuracy. Independent opening of manometer is prohibited. Repair the Device only in authorized organizations.
- Service life of device is indicated in section TECHNICAL SPECIFICATIONS. Service life is determined from the moment of commodity delivery to the customer. On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device.
- Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.
- The arm cuff may withstand multiple sanitary treatments. The internal tissue surface of the arm cuff (contacting on arm) may be cleaned with cotton ball soaked in 3%-solution of hydrogen peroxide. At long use the partial color fading of the tissue coating of the arm cuff is possible. Washing and ironing of the arm cuff are not allowed.

WARRANTY

The purchased device is warranted for 24 months from the date of sale. The warranty period for the cuff and bulb is 12 months from the date of sale. Warranty obligations shall be formalized with a warranty card when the device is sold to a customer. Faults detected during the warranty period will be repaired free of charge by a specialized service center within 21 business days from the date of device delivery to the service center. The costs for the defective device transportation are borne by Little Doctor Europe Sp. z o.o. o. (door-to-door service), except for faults that are not covered by the warranty. The party submitting the claim must properly prepare the device for shipment, protecting it from damage during transportation and enclose a completed warranty card (with the seller's stamp, date of sale and seller's signature) and proof of purchase (invoice or receipt). The warranty does not cover damages resulting from inadequate protection of the device during transportation. The warranty period is extended for the time the device is under repair. Post-warranty repair of the device is performed at the user's expense. The warranty does not cover the following: consumables and accessories for individual use (tubes, nebulizers, masks, mouthpieces, etc.); products with mechanical, thermal or chemical damage, or damage caused by violations of the rules specified in the operating instructions; products with signs of modification, opening and/or repair by an unauthorized service center (by an individual); products that have been connected to unrecommended power supplies or used with unrecommended consumables; products that have been operated with non-remedied defects; products that are defective through no fault of the manufacturer, such as: actions of third parties, natural phenomena and natural disasters, animals, insects, foreign objects or liquids entering the device, etc.; internal and external contamination, scratches, cracks, scrapes and other mechanical damage to the equipment occurring during its operation; products with illegible or altered serial numbers or a damaged warranty seal; batteries, glass parts, light bulbs, as well as operations such as adjusting, cleaning and other maintenance of the product. The buyer has the right to replace the device with a serviceable one if an authorized service center finds a manufacturing defect that cannot be eliminated or after 3 unsuccessful repairs of the same malfunction.

TECHNICAL SPECIFICATIONS	
Pressure measurement range, mmHg	from 0 to 300 (pressure in an arm cuff)
Pressure indication range, mmHg	from 0 to 300
Range of admissible absolute error at measurement of air pressure in an arm cuff, mmHg	±3
Operation conditions:	
Temperature, °C	from 10 to 40
Relative humidity, % Rh	85 and lower
Storage and transportation conditions:	
Temperature, °C	from minus 34 to 65
Relative humidity, % Rh	85 and lower
Weight (without package, case):	
model LD-70 / LD-70NR, g	270/240
model LD-71, LD-71A, g	330
model LD-60 / LD-61, g	360/310
model LD-80, g	330
model LD-100, g	480
Dimensions (package), mm:	
model LD-60, LD-61, LD-70, LD-70NR, LD-71, LD-71A, LD-80	115 x 185 x 75
model LD-100	173 x 135 x 140
Year of manufacture	The year and month of manufacture are indicated on the lower part of the manometer in the serial number after the letter "A"

CUFF SIZE OF THE DEVICE:	
model LD-70 model LD-71 model LD-71A model LD-100	adult's for upper arm circle from 25 to 36 cm
model LD-61 model LD-70NR	adult's for upper arm circle from 25 to 40 cm
model LD-60	big adult's for upper arm circle from 33 to 46 cm
model LD-80	for newborn infants: for upper arm circle from 7 to 12 cm for infants for upper arm circle from 18 to 26 cm

SYMBOL EXPLANATION:
C 0123 Compliance with the Directive MDD 93/42/EEC
i Important: Read the manual
EU REP Authorized Representative
M Manufacturer
I Importer
UDI Unique Device Identification
Protect from moisture
Storage/operation temperature and humidity
SNI Serial number
MD Medical Device
Date of manufacture

This Device manufacturing is certified according to international standard ISO 13485. Device complies with the requirements of European Directive MDD 93/42/EEC. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Revision date of the present Manual is indicated on the last page as XX-XXXX-YYMM-NN, where YY is the year, MM is the month and NN is the number of revision.

INFORMATION ON THE MANUFACTURER, AUTHORIZED REPRESENTATIVE AND IMPORTER

✉ **Complaints and requests should be addressed to:**
Little Doctor Europe Sp. z o.o. 57G Zawila Street, 30-390, Krakow, Poland
service phone: +48 12 2684748, 2684749

Manufactured under contract for Little Doctor International (S) Pte. Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE 199591. Postal Address: 3 Coleman Street #03-24 Singapore 179804).

Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No.8, Tongxing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).

EU REP: Little Doctor Europe Sp. z o.o., 57G Zawila Street Krakow 30-390 Poland.
Little Doctor Europe Sp. z o.o., 57G Zawila Street Krakow 30-390 Poland.

For more information please visit www.littledoctor.org

ADDITIONAL ACCESSORIES

Cuffs "LD-Cuff" of the following sizes may be acquired in addition to available cuffs:

Name / Size	Upper arm circle, cm	Material	Pneumatic chamber	
C2N / for newborn infants	7-12	cotton	LD-S02N	without ring
C2I / for infants	11-19	cotton	LD-S02I	
C2C / child's N2C / child's	18-26	cotton nylon	LD-S02C	
C2A / adult's N2A / adult's	25-40	cotton nylon	LD-S02A	
C2L / big adult's	34-51	cotton nylon	LD-S02L	with ring*
C2T / for upper leg	40-66	cotton	LD-S02T	
N2AR / adult's	25-36	nylon	LD-S02A	
N2LR / big adult's	33-46	nylon	LD-S02L	

* Metal fixing ring on cuff for convenient independent putting on of cuff.

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